

PST-60HL, Plate Shaker-Thermostat

Premium Product Class

DESCRIPTION

Plate Thermo-Shakers are designed for shaking and thermostating 2 standard 96-well microplates.

A multisystem principle, used in design of the Thermo-Shaker, allows operating it as 3 independent devices:

- Incubator;
- · Microplate shaker;
- Thermo-Shaker.

A distinctive feature of Biosan Plate Thermo–Shakers is thepatented by the company Two–Side Microplates Heating, which allows to achieve full correspondence of the set and actual temperature in the microplate wells.

Thermo-shaker provides heating up to 60°C, which is sufficient for carrying out ELISA tests.

Plate Shaker-Thermostat provides:

- Soft or intensive sample shaking
- Rotation speed regulation, stabilization and indication
- Even rotation amplitude throughout the Thermo-Shaker platform
- Required operation time setting and indication
- Automatic stopping of the platform movement after expiration of the set time
- Setting and indication of the required temperature on the platform
- Automatic fault diagnostics (temperature sensor, platform heating, lid heating etc.)

Application fields:

- Cytochemistry for in situ reactions
- Immunochemistry for immunofermentative reactions
- Biochemistry for enzyme and protein analysis
- Molecular biology for micro array analysis

Temperature Calibration Function

With the help of the temperature calibration function the user can calibrate the unit approx. ±6% of the selected temperature to compensate differences in the thermal behaviour of plates from different manufacturers.



CAT. NUMBER

BS-010119-AAI	230VAC 50/60Hz Euro plug
BS-010119-AAQ	230VAC 50/60Hz UK plug
BS-010119-AA4	230VAC 50/60Hz AU plug
BS-010119-AAJ	100VAC 50/60Hz US plug, 120VAC 60Hz US plug
BS-010119-BK	IQ OQ document
BS-010119-CK	PQ document

PST-60HL, Plate Shaker-Thermostat Page 1 of 2

SPECIFICATIONS

Temperature setting range	+25°C +60°C
Temperature control range	+5°C above ambient +60°C
Temperature setting resolution	0.1°C
Temperature stability	±0.1°C
Temperature uniformity at +37°C	±0.25°C
Temperature calibration coefficient range	0.9361.063 (± 0.063)
Heating	Patented two–side microplate heating
Speed control range	250-1200 rpm (increment 10 rpm)
Digital time setting	1 min–96 hrs / non–stop* (increment 1 min)
Orbit	2 mm
Display	LCD, 16 x 2 signs
Max. height of microtest plate	18 mm
Number of microtest plates	2
Platform dimensions (w x d)	250 x 150 mm
Overall dimensions (W×D×H)	270 x 260 x 125 mm
Weight	6.1 kg
Input current/power consumption	12 V DC, 3.3 A / 40 W
External power supply	Input AC 100–240 V; 50/60 Hz; Output DC 12 V
* Timer range can be reprogrammed on customers demand	+

PST-60HL, Plate Shaker-Thermostat

Page 2 of 2