

V-1 plus, Personal Vortex

DESCRIPTION

Vortex V-1 plus is an ideal instrument for gentle mixing to vigorous resuspension of cells and biological and chemical liquid components in tubes using eccentric mechanism.

Vortex has two modes:

- 1. Continuous operation;
- $2. \ \ Impulse operation (activated by pressing the cap with the tube's bottom).$

SPECIFICATIONS

Eccentric mixing principle	+
Speed control range	500-3000 RPM
Acceleration time	<1 s
Maximum continuous operation time	24 h
Mixing module for tubes	from 0.2 to 50 ml
Maximum mixing volume	30 ml
Orbit	4 mm
Overall dimensions (W×D×H)	90x150x80 mm
Weight	0.8 kg
Input current/power consumption	12 V, 320 mA / 3.8 W
External power supply	Input AC 100–240 V; 50/60 Hz; Output DC 12 V





CAT. NUMBER

BS-010203-AAG	230VAC 50/60Hz Euro plug
BS-010203-AAK	230VAC 50/60Hz UK plug, 230VAC 50/60Hz AU plug, 100VAC 50/60Hz US plug, 120VAC 60Hz US plug
BS-010203-BK	IQ OQ document
BS-010203-CK	PQ document

V-1 plus, Personal Vortex Page 1 of 1