

# JT Test after Django upgrade

## DESCRIPTION

JT Test after Django upgrade

## CAT. NUMBER

jt-100

jt 100

## SPECIFICATIONS

|                           |   |
|---------------------------|---|
| Processing capability     | 1–12 samples per run                    |
| Tube's volume             | 1.5 - 2 ml                              |
| Extraction technology     | magnetic particle separation technology |
| Tube's manufacturer       | Eppendorf or equivalent                 |
| Protocol                  | on board, preinstalled protocols        |
| Protocol input            | barcode scanner                         |
| Sample volume             | 100–1000 µl                             |
| Elution volume            | 50–300 µl                               |
| Connection to PC          | not required                            |
| Display                   | LCD (20×4)                              |
| Certification             | CE IVD                                  |
| Nominal operating voltage | 110–240 V, 50/60 Hz                     |
| Dimensions (W×D×H)        | 480×700×520 mm                          |
| Weight                    | 43 kg                                   |

## ACCESSORIES



#### ScienCell products

ScienCell Research Laboratories ([ScienCell](#)) is an expanding biotechnology company whose mission is the research and development of cell and cell-related products for experimental and therapeutic use. The ScienCell ...

[read more](#)



**Bio PP-4**  
BS-010116-AK  
platform

Flat platform with non-slip silicone mat for Petri dishes, culture flasks, agglutination.